

Maintenance kit for Waters Alliance system, with gradient LPG pump, 2487 UV detector and autosampler:

Part number: MK-WAT-ALL-U-A

2 × PA-WAT-601	- Piston assembly
2 × CV-WAT-608	- Check valve cartridge
2 × WS-WAT-715	- Wash seal - Face
4 × WS-WAT-815	- Wash seal tube kit
2 × WS-WAT-615	- Wash seal - Plunger
2 × PS-WAT-602TW	- Piston seal yellow
2 × FS-WAT-915	- Face seal
1 × WAT005139	- In line filter
1 × CTS-3109	- Syringe
1 × LD-WAT-105LL	- Deuterium lamp

Part number: MK-WAT-ALL-U-B

2 × PA-WAT-601	- Piston assembly
2 × CV-WAT-608	- Check valve cartridge
2 × WS-WAT-715	- Wash seal - Face
4 × WS-WAT-815	- Wash seal tube kit
2 × WS-WAT-615	- Wash seal - Plunger
2 × PS-WAT-602TW	- Piston seal yellow
2 × FS-WAT-915	- Face seal
1 × WAT005139	- In line filter
1 × CTS-3109	- Syringe
1 × LD-WAT-105LL	- Deuterium lamp long life
1 × CTS7550	- Autosampler Needle
1 × CTS4555	- Autosampler fluid pack
1 × CTS10839	- Seal fluid pack (pump)

Maintenance kit for Waters Alliance system, with gradient LPG pump, 2996/996 DAD detector and autosampler:

Part number: MK-WAT-ALL-D-A

2 × PA-WAT-601	- Piston assembly
2 × CV-WAT-608	- Check valve cartridge
2 × WS-WAT-715	- Wash seal - Face
4 × WS-WAT-815	- Wash seal tube kit
2 × WS-WAT-615	- Wash seal - Plunger
2 × PS-WAT-602TW	- Piston seal yellow
2 × FS-WAT-915	- Face seal
1 × WAT005139	- In line filter
1 × CTS-3109	- Syringe
1 × LD-WAT-104LL	- Deuterium lamp long life

Part number: MK-WAT-ALL-D-B

2 × PA-WAT-601	- Piston assembly
2 × CV-WAT-608	- Check valve cartridge
2 × WS-WAT-715	- Wash seal - Face
4 × WS-WAT-815	- Wash seal tube kit
2 × WS-WAT-615	- Wash seal - Plunger
2 × PS-WAT-602TW	- Piston seal yellow
2 × FS-WAT-915	- Face seal
1 × WAT005139	- In line filter
1 × CTS-3109	- Syringe
1 × LD-WAT-104LL	- Deuterium lamp long life
1 × CTS7550	- Autosampler Needle
1 × CTS4555	- Autosampler fluid pack
1 × CTS10839	- Seal fluid pack (pump)